

Gen2Epi GUI Quick Tutorial

“Filename”, ‘foldername’, menu and submeu items, Open files and confirm

Locations and name of the test data used in this tutorial:


Test data directory: *‘/home/gen2epi/Desktop/Test_DATA’*

FASTQ files: *‘/home/gen2epi/Desktop/Test_DATA/WHO_Data’*


Genome files: *‘/home/gen2epi/Desktop/Test_DATA/WHO_Full_Reference_genome/Chromosome’*

Annotation files: *‘/home/gen2epi/Desktop/Test_DATA/WHO_Genome_Annotation/Chromosome’*


- 1) Click on the Gen2EpiGUI icon on the Desktop to start the program.




- 2) To use the program please follow the instruction as given in **Introductory Demo** under **Help** menu on the top right-hand side.




- 3) Load the text file (tab-separated text file listing the names of all FASTQ files), named “*Input*”, by clicking on the **Open Project** under **Project** menu and browse (and **open**) the file as shown in the picture below. File “*Input*” is present in the “*Test_DATA*” folder at the Desktop folder.




You will see the following ****SYSTEM**** message in the text window when you click on **open** which shows that the “*Input*” file is successfully loaded.




- 4) Set the output directory by clicking on **Set Output Directory** under the **Project** menu. Browse the path (navigate to where you want to save the output files generated by “Gen2EpiGUI” e.g. in this case the path is ‘/home/gen2epi/Desktop/Test_DATA/’) and click **open**.




- 5) Write the name of your output directory (e.g. “*Output1*”) in the search box and click **confirm**.



- 6) Set input directories for FASTQ files by clicking on **Set Input Directories** -> **Fastq** under the **Project** menu. Browse the folder “*WHO_Data*” and click **open** as shown in the figure below:




- 7) Set input directories for genome files by clicking on **Set Input Directories** -> **Genome File**. Select the folder (with genome files in it) by browsing the path *'/home/gen2epi/Desktop/Test_DATA/WHO_Full_Reference_genome/Chromosome'* and click **open** as shown in the figure below:




Please Note: You will not be able to see anything when you choose a directory as shown above but that's all right.

- 8) Set input directories for annotation files by clicking on **Set Input Directories** -> **Annotation** menu. Select the folder (with annotation files in it) by browsing the path *'/home/gen2epi/Desktop/Test_DATA/WHO_Genome_Annotation/Chromosome'* and click on **open** as shown in the figure below:




Again, you will not be able to see anything under your chosen directory as shown above but that's all right.


- 9) Once all paths are set properly you will see following ****SYSTEM**** MESSAGES in “Gen2EpiGUI” main text window.




- 10) Now to clean your data: Under **Data Cleaning** menu – click on **Quality Control** submenu to check the quality of your raw fastq files.



Start of the analysis



End of the analysis: On Successful completion, you will see following message



Note: Quality control results are present under *'/home/gen2epi/Desktop/Test_DATA/Output1'*. There will be two folders *"MultiQC-Raw"* and *"QualityControl"*. Please open the .html files under these folders to visualize the quality of each sample.

- 11) In order to trim the raw reads, users have to click on the **trimming** submenu that will further ask for the confirmation of the parameters. Users have the option to use either the default parameters or change them according to their requirements. After confirmation, the program will start trimming the reads.

Gen2EpiGUI - PROJECT OPEN: "/home/gen2epi/Desktop/Test_DATA/Input"

Project	Data_Cleaning	Assembly	AMR & Strain	View Output Files
Started a	Read Mapping	Ctrl-r	fastq.gz	
Approx 5			_001.fastq.gz	
Approx 1	Quality Control	Ctrl-q	2_001.fastq.gz	
Approx 1	Trimming	Ctrl-i	2_001.fastq.gz	
Approx 2	Both	Ctrl-b	2_001.fastq.gz	
Approx 2			2_001.fastq.gz	
Approx 3	Reset Trimming Parameters		2_001.fastq.gz	
Approx 3			2_001.fastq.gz	
Approx 4			2_001.fastq.gz	

↓

Project Data_Cleaning Assembly AMR & Strain View Output Files

Please confirm parameters to be used:

Leading Length [5] Trailing Length [3] Sliding Window Start [4] Sliding Window End [15] Minimum Length [30]

Confirm Cancel

↓

Start of the analysis

Project Data_Cleaning Assembly AMR & Strain View Output Files

```

**SYSTEM** Executing Trimming script, please wait...
perl /home/gen2epi/Desktop/Gen2Epi_Scripts/WGS_SIBP_P1.pl /home/gen2epi/Desktop/Test_DATA/Input /home/gen2epi/Desktop/Test_DATA/WHO_Data Trimming /home/gen2epi/Desktop/Test_DATA/Output1 3 3 4:15 30 2>=1
TrimmomaticPE: Started with arguments: -phred33 /home/gen2epi/Desktop/Test_DATA/WHO_Data/WHO-F_S2_L001_R1_001.fastq.gz /home/gen2epi/Desktop/Test_DATA/WHO_Data/WHO-F_S2_L001_R2_001.fastq.gz /home/gen2epi/Desktop/Test_DATA/Output1/Trimming/OutputPaired_WHO-F_S2_L001_R1_001.fastq.gz /home/gen2epi/Desktop/Test_DATA/Output1/Trimming/OutputUnpaired_WHO-F_S2_L001_R1_001.fastq.gz /home/gen2epi/Desktop/Test_DATA/Output1/Trimming/OutputPaired_WHO-F_S2_L001_R2_001.fastq.gz /home/gen2epi/Desktop/Test_DATA/Output1/Trimming/OutputUnpaired_WHO-F_S2_L001_R2_001.fastq.gz LEADING:3 TRAILING:3 SLIDINGWINDOW:4:15 MINLEN:30
Multiple cores found: Using 4 threads

```

End of the analysis: On Successful completion, you will see following message


```

Searching 8 files.
[INFO] fastqc : Found 8 reports
[INFO] multiqc : Compressing plot data
[INFO] multiqc : Report : Desktop/Test_DATA/Output1/MultiQC-Trimmed/multiqc_report.html
[INFO] multiqc : Data : Desktop/Test_DATA/Output1/MultiQC-Trimmed/multiqc_data
[INFO] multiqc : MultiQC complete
Finish Data Cleaning
**SYSTEM** Script completed (Trimming)...
**SYSTEM** Logged script output to /home/gen2epi/Desktop/Test_DATA/Output1/SystemLogs/Trimming.log


```

Note: Trimming results are present under ‘/home/gen2epi/Desktop/Test_DATA/Output1’. There will be three folders “MultiQC-Trimmed”, “Trimming” and “Trimmed_QC”. Please open the .html files under these folders to visualize the quality of each sample.


- 12) It is also possible to run step 10 and 11 together as a single command. First open ‘/home/gen2epi/Desktop/Test_DATA/Output1’ folder and delete the “MultiQC-Raw”, “quality control”, “MultiQC-Trimmed”, “Trimming” and “Trimmed_QC” folders. Now go to “Gen2EpiGUI” and click on the **both** tab under **Data Cleaning** menu.



Start of the analysis




End of the analysis: On Successful completion, you will see following message



Note: - You will find all the results under '/home/gen2epi/Desktop/Test_DATA/Output1' in five different folders i.e. "MultiQC-Raw", "quality control", "MultiQC-Trimmed", "Trimming" and "Trimmed_QC"

13) Now, perform the *de novo* assembly of trimmed reads by clicking on **De Novo** under **Assembly** menu. You will see the output as shown in the picture below.

Please Note: - This step may take a while depending on individual computer configurations. It could be possible that nothing happens when you try to press enter in the VM image, no need to worry, VM will resume itself after completing the *de novo* assembly part. You can minimize the VM image window and come back later.



14) Once the *de novo* assembly (step 13) is completed. You will see the following message

```
===== Mismatch correction finished.
* Corrected reads are in /home/gen2epi/Desktop/Test_DATA/Output1/Plasmid_AssemblyTrimmedReads/WHO-L/corrected/
* Assembled contigs are in /home/gen2epi/Desktop/Test_DATA/Output1/Plasmid_AssemblyTrimmedReads/WHO-L/contigs.fasta
* Assembled scaffolds are in /home/gen2epi/Desktop/Test_DATA/Output1/Plasmid_AssemblyTrimmedReads/WHO-L/scaffolds.fasta
* Assembly graph is in /home/gen2epi/Desktop/Test_DATA/Output1/Plasmid_AssemblyTrimmedReads/WHO-L/assembly_graph.fastg
* Assembly graph in GFA format is in /home/gen2epi/Desktop/Test_DATA/Output1/Plasmid_AssemblyTrimmedReads/WHO-L/assembly_graph_with_scaffolds.gfa
* Paths in the assembly graph corresponding to the contigs are in /home/gen2epi/Desktop/Test_DATA/Output1/Plasmid_AssemblyTrimmedReads/WHO-L/contigs.paths
* Paths in the assembly graph corresponding to the scaffolds are in /home/gen2epi/Desktop/Test_DATA/Output1/Plasmid_AssemblyTrimmedReads/WHO-L/scaffolds.paths

===== SPAdes pipeline finished.

SPAdes log can be found here: /home/gen2epi/Desktop/Test_DATA/Output1/Plasmid_AssemblyTrimmedReads/WHO-L/spades.log


Thank you for using SPAdes!

**SYSTEM** Script completed (De_Novo)...


**SYSTEM** Logged script output to /home/gen2epi/Desktop/Test_DATA/Output1/SystemLogs/De_Novo.log
```

Note: You will find the results under ‘/home/gen2epi/Desktop/Test_DATA/Output1’ in four folder i.e. “Chrom_AssemblyTrimmedReads” (assembled chromosome contigs in FASTA format generated from trimmed reads), “ChromContigAssemblyTrimmedStat” (chromosome assembly statistics), “Plasmid_AssemblyTrimmedReads” (assembled plasmid contigs in FASTA format generated from trimmed reads), and “PlasmidContigAssemblyTrimmedStat” (plasmid assembly statistics)

15) Perform the scaffolding of *de novo* assembly by clicking on **Assembly -> Scaffolding -> Chromosome -> Multiple Genomes** as shown below:



Start of the Scaffolding process



```
Gen2EpiGUI - PROJECT OPEN: "/home/gen2epi/Desktop/Test_DATA/Input"
Project Data Cleaning Assembly AMR & Strain View Output Files
**SYSTEM** Executing Multiple_Genomes script, please wait...
perl /home/gen2epi/Desktop/Gen2Epi_Scripts/WGS_SIBP_P3-Chr-C1.pl /home/gen2epi/Desktop/Test_DATA/Input/home/gen2epi/Desktop/Test_DATA/WHO_Full_Reference_genome/Chromosome /home/gen2epi/Desktop/Test_DATA/Output1/Chrom_AssemblyTrimmedReads /home/gen2epi/Desktop/Test_DATA/WHO_Genome_Annotation/Chromosome 2 TXT /home/gen2epi/Desktop/Test_DATA/Output1 2-&1
[09:35:34] INFO: Starting Ragout v2.0
[09:35:34] INFO: Synteny block scale set to 'small'
[09:35:34] INFO: Running Sibelia with block size 5000
```


Scaffolding process completion




```
Finished: 2019-06-03 09:41:13
Elapsed time: 0:00:07.934315
Total NOTICES: 5; WARNINGS: 5; non-fatal ERRORS: 0
Thank you for using QUAST!
Attaching package: &#x26;#x26;dplyr&#x26;#x26;
The following objects are masked from &#x26;#x26;package:stats&#x26;#x26;:
  filter, lag
The following objects are masked from &#x26;#x26;package:base&#x26;#x26;:
  intersect, setdiff, setequal, union
Joining, by = "samples"
**SYSTEM** Script completed (Multiple_Genomes)...
**SYSTEM** Logged script output to /home/gen2epi/Desktop/Test_DATA/Output1/SystemLogs/Multiple_Genomes.log
```

Note: Output files generated from this step will be stored in folder “*Chr_Scaffolds*” under ‘/home/gen2epi/Desktop/Test_DATA/Output1’.

- 16) To identify the different types of *Neisseria gonorrhoeae* plasmids using assembled plasmid contigs (in step 13), select **Assembly -> Scaffolding -> Plasmid**.



Once you click on **Plasmid**. You will see the following output




```
Gen2EpiGUI - PROJECT OPEN: "/home/gen2epi/Desktop/Test_DATA/Input"
Project Data Cleaning Assembly AMR & Strain View Output Files
**SYSTEM** Executing Scaffolding_Plasmid script, please wait...
perl /home/gen2epi/Desktop/Gen2Epi_Scripts/WGS_SIBP_P3-Plas_C1.pl /home/gen2epi/Desktop/Test_DATA/Input/home/gen2epi/Desktop/Test_DATA/Output1/Plasmid_AssemblyTrimmedReads 1 /home/gen2epi/Desktop/Test_DATA/Plasmid.fasta /home/gen2epi/Desktop/Test_DATA/Output1 2-&1
**SYSTEM** Script completed (Scaffolding_Plasmid)...
**SYSTEM** Logged script output to /home/gen2epi/Desktop/Test_DATA/Output1/SystemLogs/Scaffolding_Plasmid.log
```


Note: Output generated from this step will be stored in folder “*Plasmid_Identification*” under ‘/home/gen2epi/Desktop/Test_DATA/Output1’.

17) NgMAST Typing

- a. Update the underlying MAST database by clicking on **AMR & Strain -> MAST -> MAST Update** submenu.




MAST database update




- b. Perform the MAST typing by selecting **AMR & Strain -> MAST -> MAST Typing** submenu



Strain typing




Note: MAST strain typing for each sample is present in file "NgMAST.txt" under '/home/gen2epi/Desktop/Test_DATA/Output1' as shown below.




18) NgMLST Typing


- Update the underlying MLST database by clicking on **AMR & Strain -> MLST -> MLST Update** submenu.



MLST database update



- Perform the MLST typing by selecting **AMR & Strain -> MLST -> MLST Typing** submenu.



Start of the strain typing analysis

```
**SYSTEM** Executing MLST_Typing script, please wait...
perl /home/gen2epi/Desktop/Gen2Epi_Scripts/WGS_SIBP_P4_Epi.pl /home/gen2epi/Desktop/Test_DATA/Input /home/gen2epi/Desktop/Test_DATA/Output1/Chr_Scaffolds MLST /home/gen2epi/Desktop/Test_DATA/Output1 /home/gen2epi/Desktop/Gen2Epi_Scripts/MLST-Genes.fasta /home/gen2epi/Desktop/Gen2Epi_Scripts/MLST_alleles.fasta /home/gen2epi/Desktop/Gen2Epi_Scripts/pubMLST_profile.txt 2>&1
```

Completion of the strain typing analysis

```
Building a new DB, current time: 06/03/2019 10:38:30
New DB name: /home/gen2epi/Desktop/Gen2Epi_Scripts/MLST_alleles.fasta
New DB title: /home/gen2epi/Desktop/Gen2Epi_Scripts/MLST_alleles.fasta
Sequence type: Nucleotide
Deleted existing Nucleotide BLAST database named /home/gen2epi/Desktop/Gen2Epi_Scripts/MLST_alleles.fasta
Keep MBits: T
Maximum file size: 1000000000B
Adding sequences from FASTA; added 7219 sequences in 0.144738 seconds.

**SYSTEM** Script completed (MLST_Typing)...

**SYSTEM** Logged script output to /home/gen2epi/Desktop/Test_DATA/Output1/SystemLogs/MLST_Typing.log
```

Note: MLST strain typing for each sample is present in file “NgMLST.txt” under ‘/home/gen2epi/Desktop/Test_DATA/Output1’ as shown below.

Chrom_ AssemblyTrimmedReads ChromContigAssemblyTrimmedStat Chr_Scaffolds GenomeStateAll.txt MultiQC-Raw

MultiQC-Trimmed NgMAST.txt **NgMLST.txt** Plasmid_ AssemblyTrimmedReads PlasmidContigAssemblyTrimmedStat


Open NgMLST.txt ~/Desktop/Test_DATA/Output1 Save

Sample	ST	abcZ	adk	aroE	fumC	gdh	pdhC	pgm	clonal_complex
WHO-F_scaffolds_pubMLST.txt				10934	200	39	67	157	148 153 65
WHO-G_scaffolds_pubMLST.txt				1903	126	39	67	157	148 153 65
WHO-K_scaffolds_pubMLST.txt				7363	59	39	67	78	148 153 65
WHO-L_scaffolds_pubMLST.txt				1590	126	39	67	78	149 153 65


Plain Text Tab Width: 8 Ln 1, Col 1 INS

19) NgSTAR


- a. Update the underlying STAR database by clicking on **AMR & Strain** -> **STAR** -> **STAR Update** -> **STARdbUpdate** submenu.




STAR database update




- b. Update the underlying STAR database metadata by clicking on **AMR & Strain** -> **STAR** -> **STAR Update** -> **STARdbMetaDataUpdate** submenu.




STAR database metadata update




- c. Perform the STAR typing by selecting **AMR & Strain** -> **STAR** -> **STAR Typing** submenu.




Start of the STAR typing process



Completion of the STAR typing process



Note: NgSTAR output for each sample is present in file “*NgStarSearchResults-WithoutST.txt*” under ‘/home/gen2epi/Desktop/Test_DATA/Output1’ as shown below.




20) Tetracycline Resistance


- 1) Click on TETRES -> Nucleotide submenu under AMR & Strain menu to extract the rpsJ nucleotide sequences. These fasta sequences can be visualized under any multiple sequence alignment program.




Starting




Completion




Note: You can see the rpsJ nucleotide sequence for each sample in file “*Nucl_rpsJ.fasta*” under ‘/home/gen2epi/Desktop/Test_DATA/Output1’ as shown below.





- 2) Click on **TETRES** -> **protein** submenu under **AMR & Strain** menu to extract the rpsJ protein alignment. These fasta sequences can be visualized under any multiple sequence alignment program.



Note: You can see the rpsJ nucleotide sequence for each sample in the file “*Prot_rpsJ.fasta*” under ‘/home/gen2epi/Desktop/Test_DATA/Output1’ as shown below.



21) View Output Files: - Finally, the text and HTML outputs generated from the above steps are accessible under **View Output File** tab.



Gen2EpiGUI

Project Data Cleaning Assembly AMR & Strain View Output Files

Samples	NgMAST POR	TBPB		
WHO-F_scaffolds.fasta	3303	2002	165	
WHO-G_scaffolds.fasta	621	90	18	
WHO-K_scaffolds.fasta	1424	917	10	
WHO-L_scaffolds.fasta	1422	915	21	

Visualize the results